25th INTERNATIONAL BIOLOGY OLYMPIAD

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
Note! All the DNA electrophoresis instruments will be turned on simultaneously by the assistants **50 minutes after the start of the test.**

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Make sure that you have placed your samples in the agarose gel in accordance to the instructions on the Question paper before the power supply is turned on. After this, you will **NOT** be allowed to run the gel. To remind you, the assistant will ring the bell three times. Bell 1 will remind you that the electrophoresis will be carried out in 20 minutes. Bell 2 is a warning that electrophoresis will begin in two minutes' time. Bell 3 will mark the end of electrophoresis. At that time you will be asked to remove your gel and place it into the box provided for assessment. 파워 서플라이가 켜지기 전에, 문제지에 지시된 순서대로 모든 샘플이 젤에 로딩되었는지 확인하시오. 일단 전기영동이 시작 된 후에는 다시 시작 할 수 없다. 다시 한 번 강조하는데,

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Note : Use given reagents properly! No additional reagents will be provided.

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Part A (40 points)

Identification of Plasmids by Restriction Enzyme Analysis

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Introduction (도입)

A scientist was surprised to find unlabelled tubes in the freezer. These tubes, containing Plasmids X, Y and Z, were arbitrarily labelled Plasmid 1, Plasmid 2 and Plasmid 3. These plasmids are indistinguishable by DNA electrophoresis because they are all 3750 bp long. Since the three plasmids differ in their restriction pattern when using *EcoRI* and/or *HindIII* (Figure 1), you will conduct restriction enzyme analysis to correctly identify the three plasmids.

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
Figure 1. Restriction maps of Plasmids X, Y, Z.

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Table I in the **Answer Sheet** shows the experimental design you will perform. The experiment is divided into Series 1 and Series 2, carried out simultaneously. In each series, all three plasmids should be cut with the same restriction enzyme(s). The final concentration of the restriction buffer solution should be 1X.

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Question 1.1 Determine how these restriction enzyme(s) can be used to differentiate between the plasmids. Complete **Table I** on the **Answer Sheet** to indicate your chosen protocol for each series! (8 points). You must choose one or two enzymes for each series that are most informative based on Figure 1. If the enzyme is used in the reaction, add 1 µL enzyme.

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Question 1.2. Complete the figure on your **Answer Sheet** to show the expected molecular weight(s) of the completely digested plasmid for Plasmids X, Y and Z that you would observe using each of your two chosen Series. Indicate the restriction enzyme(s) used for Series 1 and Series 2 of your experimental design. (12 points – 2 points each lane).

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Plasmid restriction and DNA electrophoresis protocol

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Note: Electrophoresis will start 50 minutes after the test begins. After this, you will NOT be allowed to run the gel.

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Note: Spin down all reagents in microtubes before pipetting (Be sure to balance the minicentrifuge by placing the microtubes opposite each other, even if there is only one tube to spin).

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1. Label the microtubes S1 to S6, and prepare the reaction mixtures according to Table I in the Answer Sheet. Spin down the mixture by placing the microtubes in the minicentrifuge for 10-20 seconds.

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3. Add **ONLY 2 μL** of DNA staining solution into microtube labelled 1 Kb Ladder. Spin down the solution for 10 seconds.

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4. After incubation, raise your bench number and your tubes will be returned to you. Add **ONLY 1 μL** of DNA staining solution into each of these tubes, mix well by pipetting and spin down any residual liquid.

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5. Open the gel agarose package and place the gel with its tray in the electrophoresis chamber. Pour the entire TAE buffer from the bottle into the tank.

Note: Black plug is cathode (-) and red plug is anode (+).

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6. Load 10 μl of samples S1 to S6 and molecular weight marker solution into the agarose gel as shown in **Figure 2**.

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Figure 2

7. Close the lid of the gel chamber and notify the assistant by raising your hand that you are ready. The assistant will connect the cables to the power supply.

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During the electrophoresis continue with Part B.

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8. After the electrophoresis (indicated by Bell 3 ringing), the assistants will turn off the power supply. Carefully place the gel with the tray in the plastic box provided. **Pour the running buffer from the electrophoresis chamber into the plastic box.** Close the box, label the sticker with your **Student ID** and affix it to the side of the box. Leave the box on your bench. Later the gel will be documented and the picture will be attached onto your answer sheet by the assistants.

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Part B (24.5 points)
Cell Reproduction and Telomere Analysis of *Paramecium*

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Two of the phenomena displayed by *Paramecium tetraurelia* are binary fission (asexual reproduction) and conjugation (exchange of genetic materials), which are affected by the abundance of nutrients in the media. Figure 3 shows microscopic observations of two *P. tetraurelia* cultures grown in rich and poor media, respectively.

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Question 2.1. (2 points) Which of the phenomena are shown by the two cultures in Figure 3A and 3B?

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
Figure 3. *P. tetraurelia* growing in rich (A) and poor (B) media.

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Paramecium, which has just undergone conjugation, continues to asexually reproduce. In this instance, telomere shortening occurs for a certain number of divisions during logarithmic growth. **Table B-1** shows data for the growth of a strain of *Paramecium* in three replicates for four days. The cultures can be grown continuously with no significant change in the fission rate.

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Table B-1. *P. tetraurelia* growth in the first four days.

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Question 2.2. (2.5 + 8 points) Calculate the average cell concentrations and using the logarithmic value draw the growth curve of *P. tetraurelia* from Days 0 to 4 on the graph provided in your **Answer Sheet** (Precision: integer).

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Figure 4 shows the result of Southern blot analysis of *P. tetraurelia* telomeres for a total of 30 synchronized generations of cultivation. DNA from *P. tetraurelia* was digested at the beginning of the telomere region, blotted onto a membrane and then hybridized with telomere probes. Telomeres were seen as smears, the midpoints of which can be matched with molecular weight markers to determine the average telomere length.

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
Figure 4. Southern blot of *P. tetraurelia* telomeres. Lanes 1 to 5 show telomeres after 4, 7, 17, 23, and 30 generations, respectively. DNA size markers (M) are shown next to Lane 5 (in kbp).

Internal: Internal part of the chromosomal DNA that contains telomeric sequence

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Question 2.3. (12 points) Indicate in the answer sheet with a tick (✓) if the following statements are **true** or **false**.

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