Country: _____

Student Code: _____

23rd INTERNATIONAL BIOLOGY OLYMPIAD

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Introduction

Genetic mapping is routinely used in analysing the order and the identities of DNA fragments. This technique is based on the unique profiles of DNA fragments generated after DNA digests with specific combination of restriction endonucleases (RE) and revealed by DNA gel electrophoresis. It is extremely powerful for gene cloning, studying gene function and regulation, for finding candidate genes for diseases and their diagnosis and also as a forensic tool.

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Part A. Confirmation of insertion of human DNA in a cloning plasmid. (80 points)

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Using this technique, you are now tasked to confirm that a fragment of human DNA “X” (approximate size: 760 base pairs) has been inserted into a cloning plasmid or vector “V” (circular and approximate size: 2570 base pairs). You are required to design and carry out DNA digests by incubating DNA “T” with the restriction endonucleases (RE) by following the general protocol of incubation and electrophoresis given (details described below). After the gel electrophoresis, your results will be revealed by DNA staining (this will be performed by lab technicians), analysed and data interpreted.

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Protocol and Procedures

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1. Design your DNA digests (you may do a maximum of 4 tubes) in a total volume of 20 μl by using the Table in **the Answer Sheet**.

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Q1.1 (20 points) Record the desired amounts of reagents in your plan. One example is already given for Tube 2 in the table provided. All units are in μl .

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2. Prepare the mixtures by carefully pipetting the correct amount of the reagents and gently mix them by pipetting them up and down in each tube. Label the tubes. Do not contaminate one sample with another when preparing the mixture. Use a clean pipette tip for each operation. Note: use p10 micropipetter (white-coded) for pipetting reagents of less than 10 μl . [NOTE: there will be a penalty of 20 points if additional samples are requested. Please prepare the samples carefully.

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3. Spin down the mixture by placing all four tubes in the micro-centrifuge (please balance the spin by placing tubes opposite to each other). During preparation and after spinning, always keep the tubes on ice.

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4. After all the tubes have been prepared, remove them from the ice and place them into the labelled floatation rack and incubate them for 20 minutes (stopwatch is provided) at 37 $^{\circ}\text{C}$ in the

water bath assigned to you. Make sure that you retrieve your own samples after the 20 minutes' incubation time.

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Q1.2 (2 points × 5 = 10 points) Indicate true statement(s) with a tick (✓) and false statement(s) with a cross (✗).

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Q1.3 (2 points × 5 = 10 points) Which of the following principles is true of separating DNA by gel electrophoresis? Indicate true statement(s) with a tick (✓) and false statement(s) with a cross (✗).

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- 6. When the 20 minutes of DNA digests duration is up, retrieve your own tubes from the water bath.

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- 7. Add 4 μ l of DNA loading-dye (blue colour). Mix them well by pipetting the mixture up and down and spin down any residual liquid using the micro-centrifuge.

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- 8. Using the p100 micropipetter (yellow-coded), load 15 μ l of the sample mixture of DNA digests with the loading dye into the “wells” of the agarose gel provided. Make sure that you position the pipette tips carefully on top of the wells and gently deliver the mixture to the wells without spilling them. Load 15 μ l of each of the Markers, L1 and L2. Add your samples according to the following scheme of lanes, starting from the left end of the gel.

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- 9. Cover a gel with a lid and connect the power supply to run at 100 volts for 20 min. Please be careful and do not touch any part of the electrodes and power supply.

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
Q1.4 (20 points) Consider the following scenario: A piece of linear human DNA (1 kbp) was digested by a particular enzyme RE3, resulting in 2 fragments of 650 bp and 350 bp. The same piece of 1 kbp DNA was digested with another enzyme RE4, releasing 2 fragments of 800 bp and 200 bp. And when this 1 kbp DNA was digested with RE3 and RE4 together, 3 fragments of DNA were generated, 650 bp, 200 bp and 150 bp.

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Sketch a linear map of this piece of DNA by indicating the position of RE3 and RE4 digests in the space provided. An example of such a sketch is provided below as a guide.

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12. When the 20 minutes of gel running time is up, turn off the power supply and remove the lid of the gel tank. Carefully remove the gel (still on the gel tray) and place it on the petri dish provided. Bring your gel to the gel doc that has been assigned for your usage and the technician will photograph it for you.

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Q1.5 (10 points) Your skills in running a gel will be assessed by the quality of the gel produced.

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
Q1.6 (1 point × 5 = 5 points) Using the DNA ladder markers (in basepairs) provided below as the reference, estimate the sizes of the fragments/bands. You may draw a line across the band of your query and the size marker to do the estimation. How many fragment(s) of DNA were generated by RE1 and RE2? And what is/are the estimated size(s)? Answer using numerals.

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Q1.7 (1 point) What is the estimated size of the test DNA sample (T)? Answer using numerals.

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Q1.8 (1 point) Based on your results, is the test DNA sample (T) larger, smaller or the same size as the empty vector? Indicate your answer with a tick (✓) in the correct box.

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Q1.9 (1 point) Does the test DNA sample (T) contain any insert? Indicate yes with a tick (✓) and no with a cross (✗).

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Q1.10 (2 points) Uncut DNA appears to move faster than any of the samples digested with RE2. Why? Indicate your answer with a tick (✓) in the correct box.

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Part B. Determination of orientation by which fragment was inserted. (20 points)

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Q1.11 (20 points) Construct possible the restriction map(s) for the DNA "T" by indicating the relative position of RE1 and RE2 and the distance between them in the Answer Sheet. An example of such a map is provided below as a guide.

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